




ADVANCING DENTAL INSTRUMENT COMPANY

IMPLANT REMOVER KIT(S)

No more stress when removing the fixture!

Only one solution

→ Tilt →	→ Accurate drillings
→ Excessive bone deletion →	→ Minimal bone deletion
→ Risk of fracture →	→ Resolving risk of fracture
→ Increasing chair time →	→ Reducing chair time



- ✓ Easy to remove the fail implant through accurate engage.
- ✓ Accurate drilling and chair time reduction.
- ✓ Minimal bone removal.
- ✓ Eliminating the risk of fracture.

EASY REMOVER

SIZE	INNER -OUTER (mm)	IMPLANT REMOVING SIZE
Ø 3.5	3.9 - 4.9	Available in sizes less than 3.9mm
Ø 4.0	4.3 - 5.3	Available in sizes less than 4.3mm
Ø 4.5	4.7 - 5.7	Available in sizes less than 4.7mm
Ø 5.0	5.2 - 6.2	Available in sizes less than 5.2mm
Ø 6.0	6.1 - 7.1	Available in sizes less than 6.1mm